



For high efficiency transformation of T&A Cloning Vector in 1-10 minutes, utilise unique RBC HIT Non-Heat Shock Competent Cells!

- HIT Competent Cells™-DH5a: Value 10⁸, High 10⁸ and Super 10⁹.
- HIT Competent Cells™-JM109: Value 10⁸ and High 10⁸.
- HIT Competent Cells™-Blue: Value 10⁸, High 10⁸ and Super 10⁹.
- HIT Competent Cells™-BL21: Value 10⁷.



Real Biotech Corporation

Copyright © 2009 Real Biotech Corporation. All Rights Reserved.
Real Biotech Corporation is a leading R&D based life science solutions provider.
For more information on our extensive and innovative life science range, please visit our website.

RBC Cloning Systems RBC T&A Cloning Kit Protocol Book

Rapid and Economical Cloning PCR Product

RC001 RBC T&A Cloning Kit / RC013 RBC T&A Cloning Vector

Contents

- **Product Description** P1
- Introduction P1
- Contents P1
- Shipping Conditions P2
- Storage Conditions P2
- Application Table P3
- Vector Map of RBC T&A Cloning Vector P4
- DNA Sequence of RBC T&A Cloning Vector Multiple Cloning Site P6
- Restriction Enzymes and RBC T&A Cloning Vector P8
- Full Length of RBC T&A Cloning Vector P12
- **Protocol for Ligation Using RBC T&A Cloning Vector** P16
- **Notes** P18
- **Troubleshooting** P20

Product Use Limitations: All products provided by Real Biotech Corporation are developed, designed and sold for research purposes only. They are not to be used for human diagnostic or drug purpose.

Product Description

Introduction

Real Biotech Corporation T&A Cloning System is ideal for rapid cloning PCR product generated using a thermostable DNA polymerase, such as RealTaq Polymerase, which adds a single terminal 3'-dA nucleotide overhang.

Contents

- Cat. No. RC001 RBC T&A Cloning Kit (20 reactions)
 - * T&A Cloning Vector (25 ng/μl): 40 μl.
 - * Control Insert DNA (10 ng/μl): 10 μl.
 - * T4 DNA Ligase (3 U/μl): 20 μl.
 - * 10X Ligation Buffer A: 100 μl.
 - * 10X Ligation Buffer B: 100 μl.
 - * Forward Primer (M13-F)(10 μM): 50 μl.
 - * Reverse Primer (M13-R)(10 μM): 50 μl.
 - * Protocol Book.

RBC T&A Cloning Kit

— Cat. No. RC013 RBC T&A Cloning Vector (20 reactions)

* T&A Cloning Vector (25 ng/μl): 40 μl.

* Protocol Book.

Shipping Conditions

RBC T&A Cloning Kits have been monitored to ensure -20°C AT ALL TIMES during batch processing, shipping and storage. RBC T&A Cloning Kits shall be shipped with dry ice.

Storage Conditions

RBC T&A Cloning Kits should be stored immediately upon receipt at -20°C in a constant temperature freezer. RBC T&A Cloning Kits can be stored for up to 12 months without showing any deduction in performance and quality with proper storage.

2

Application Table

	RBC T&A Cloning Kit	RBC T&A Cloning Vector
PCR amplified DNA (Taq)	Yes	Yes
Terminal dA Tailed (Taq)	Yes	Yes
PCR amplified DNA (Pfu)	No	No
Blunt end-dephosphorylated	No	No
Sticky end	No	No
Linker Ligation	No	No


Note:

The use of the PCR process is covered by patents obtained by Hoffman-La Roche. PCR is a trademark of Hoffman-La Roche.

3

Vector Map of RBC T&A Cloning Vector

Multiple Cloning Site	434	to	490
lacZ gene	511	to	149
AP ^r gene	2528	to	1671
T7 promoter	402	to	439
M13 universal primer	359	to	375
M13 reverse primer	528	to	507



* Before the insert incorporated into the RBC T&A Cloning Vector, there is only 1 HindIII site and no BglII site. When the incorporation is done, the T and A nucleotides on the insert will complement the sequence on the vector and generate two new sites.

4

5

DNA Sequence of RBC T&A Cloning Vector Multiple Cloning Site

The diagram shows the DNA sequence of the RBC T&A Cloning Vector Multiple Cloning Site. It includes the following elements:

- 307**: 5'-----AGCTGG CGAAAG GGGGAT GTGCTG CAAGGC GATTAA GTTGGG TAACGC CAGG
3'-----TCGACC GCTTTC CCCCTA CACGAC GTTCCG CTAATT CAACCC ATTGCG GTCC
- 359**: M13 universal primer (GTTTTCCC AGTCAC GACGTT GTAAAA CGACGG CCAGTG AATTGT AATACG ACTCAC TATAGG) and T7 promoter (CAAAGGG TCAGTG CTGCAA CATTTC GCTGCC GGTCAC TTAACA TTATGC TGAGTG ATATCC).
- 421**: Multiple cloning site with restriction enzyme sites: KpnI, SmaI, EcoRI, [HindIII], [BglII], BamHI, XbaI, Sall. The sequence is: GCGAGC TCGGTA CCCGGG CGAATTCC AAGCTT -AGATCT GGATCC CCTCTA GAGTCG CGCTCG AGCCAT GGGCCC GCTTAAGG TTCGAA- TCTAGA CCTAGG GGAGAT CTCAGC.
- 475**: ACCTGC AGGCAT GCAAGC TTGGCG TAATCA TGGTCA TAGCTG TTCCTT GTGTGA AATTGT TATCCG TGGACG TCCGTA CGTTTCG AACCGC ATTAGT ACCAGT ATCGAC AAAGGA CACACT TTAACA ATAGGC
- 541**: CTCACA ATTCCA CACAAC ATACGA GCCGGA AGCATA AAGTGT AAAGCC TGGGGT -----3'
GAGTGT TAAGGT GTGTTG TATGCT CGGCCT TCGTAT TTCACA TTCGCG ACCCCA -----5'

An "Insert DNA" box is located below the multiple cloning site.

Restriction Enzymes and RBC T&A Cloning Vector

1. No restriction sites to cut RBC T&A Cloning Vector for the following enzymes:

AatI	BbsI	BsiWI	BstBI	Eco47III	MroI	NsiI	Ppu10I	SpeI
AccIII	BclI	BsmFI	BstEII	EcoNI	MscI	NspV	PpuMI	SspBI
AflII	BfrI	BsmI	Bsu36I	EcoRV	MunI	PacI	RsrII	StyI
AgeI	BlnI	Bsp120I	CelII	EspI	NaeI	PaeR7I	SacII	Swal
ApaI	Bpu1102I	BspDI	Clal	HpaI	NcoI	PflMI	SexAI	Tth111I
AspI	BpuAI	BspEI	Csp45I	KspI	NgoMI	PinAI	SfiI	Van91I
AsuII	BsaAI	BsrGI	DrallI	MfeI	NheI	PmaCI	SfuI	XcmI
AvrII	BseAI	BssHII	EagI	MluI	NotI	PmeI	SgrAI	XhoI
BbrPI	BsgI	Bst1107I	EcXI	MluNI	NruI	PmlI	SnaBI	XmaIII

2. Single restriction site to cut RBC T&A Cloning Vector for the following enzymes:

Name	Position	Name	Position	Name	Position	Name	Position	Name	Position
AatII	2664	AspEI	1742	Cfr10I	1822	MamI	457	SspI	2546
Acc65I	430	AvaI	434	Drall	2718	NarI	237	XbaI	466
AccI	473	BanII	428	Eam1105I	1742	NdeI	185	XmaI	434
AcsI	441	BamHI	458	Ecl136II	426	PstI	482	XmnI	2341
AflIII	849	BcgI	2281	EcoO109I	2718	SacI	428		
AhdI	1742	Bpml	1812	EcoRI	441	Sall	472		
AlwNI	1265	BsaBI	457	HincII	474	SapI	733		
ApoI	441	BsaI	1803	HindII	474	Scal	2222		
Asp700	2341	BspMI	485	KasI	236	SmaI	436		
Asp718	430	BsrFI	1822	KpnI	434	SphI	488		

3. More than one restriction site to cut RBC T&A Cloning Vector for the following enzymes:

Name	Position	Name	Position	Name	Position
Acil	34sites	Cfol	17 sites	Mbol	16 sites
Alul	17sites	Csp61	169,431,2221	Mboll	7 sites
Alw44I	178,1163,2409	Ddel	6 sites	McrI	5 sites
AlwI	10sites	Dpnl	16 sites	MnII	13 sites
Aosl	259,1964	DpnII	16 sites	Msel	13 sites
ApaLI	178,1163,2409	Dral	1608,1627,2319	MslI	1994,2153,2512
Asel	620,679,1914	DrdI	98,957	MspA1I	6 sites
AsnI	620,679,1914	DsaV	12 sites	Mspl	13 sites
AspHI	5 sites	EaeI	389,688,2130	Mval	5 sites
Avall	1880,2102	EarI	297,733,2537	Mvnl	10 sites
Avill	259,1964	Eco57I	1397,2409	Mwol	13 sites
BanI	4 sites	EcoRII	5 sites	NciI	7 sites
BbvI	12 sites	Esp3I	4,46	NdelI	16 sites
Bfal	4 sites	Fnu4HI	19 sites	NlaIII	11 sites
BglI	252,1862	FnuDII	10 sites	NlaIV	11 sites
		FokI	5 sites	NspI	42,488,853

Name	Position	Name	Position	Name	Position
BmyI	5 sites	FspI	259,1964	PleI	5 sites
BsaHI	237,2279,2661	Haell	240,727,1097	Psp1406I	1968,2341
BsaJI	4 sites	HaellI	11 sites	PvuI	280,2112
BsaWI	1055,1202,2033	Hgal	4 sites	PvuII	309,673
BsiEI	5 sites	HgiAI	5 sites	RcaI	1569,2577,2682
BsiHKA1	5 sites	HhaI	17sites	RsaI	170,432,222
BsiYI	6 sites	HindIII	(448),490	Sau3AI	16 sites
BsII	6 sites	Hinfl	7 sites	Sau96I	6 sites
BsmAI	4 sites	HinPII	17 sites	ScrFI	12 sites
Bsp1286I	5sites	HpaII	13 sites	SfaNI	8 sites
BspHI	1569,2577,2682	HphI	7 sites	Sfcl	5 sites
BspWI	13 sites	HPY188I	10 sites	Snol	178,1163,2409
BsrBI	541,782,2583	ItaI	19 sites	TaqI	473,949,2393
BsrDI	1803,1977	Ksp632I	297,733,2537	TflI	684,824
Bsrl	11 sites	MaeI	4 sites	Thal	10 sites
BstNI	5 sites	Maell	5 sites	Tru9I	13 sites
BstUI	10 sites	MaellI	11 sites	Tsp509I	8 sites
BstYI	7 sites			XhoI	7 sites

Full Length of RBC T&A Cloning Vector

1	TCGGCG	GTTTCG	GTGATG	ACGGTG	AAAACC	TCTGAC	ACATGC	AGCTCC	CGGAGA	CGGTCA
	AGCGCG	CAAAAG	CACTAC	TGCCAC	TTTTGG	AGACTG	TGTACG	TCGAGG	GCCTCT	GCCAGT
61	CAGCTT	GTCTGT	AAGCGG	ATGCCG	GGAGCA	GACAAG	CCCCTC	AGGGCG	CGTCAG	CGGGTG
	ATCGAA	CAGACA	TTCGCC	TACGGC	CCTCGT	CTGTTC	GGGGAG	GGGGCG	GCAGTC	GCCCCA
121	TTGGCG	GGTGTC	GGGGCT	GGCTTA	ACTATG	CGGCAT	CAGAGC	AGATTG	TACTGA	GAGTGC
	AACCGC	CCACAG	CCCCGA	CCGAAT	TGATAC	GCCGTA	GTCTCG	TCTAAC	ATGACT	CTCACC
181	ACCATG	TGCGGT	GTGAAA	TACCGC	ACAGAT	GCGTAA	GGAGAA	AATACC	GCATCA	GGCGCC
	TGGTAT	ACGCCA	CACTTT	ATGGCG	TGTCTA	CGCATT	CCTCTT	CGTAGT	CCGCGG	
241	ATTTCG	CATTCA	GGCTGC	GCAACT	GTGGGG	AAGGGC	GATCGG	TGCGGG	CCTCTT	CGCTAT
	TAAGCG	GTAAGT	CCGACG	CGTTGA	CAACCC	TTCCCG	CTAGCC	ACGCC	GGAGAA	GCGATA
301	TACGCC	AGCTGG	CGAAAG	GGGGAT	GTGCTG	CAAGGC	GATTAA	GTGGGG	TAACGC	CAGGGT
	ATGGCG	TCGACC	GCTTTC	CCCCTA	CACGAC	GTTCCG	CTAATT	CAACCC	ATTGCG	GTCCCA
361	TTTCCC	AGTCAC	GACGTT	GTAATA	CGACGG	CCAGTG	AATTGT	AATACG	ACTCAC	TATAGG
	AAAGGG	TCAGTG	CTGCAA	CATTTT	GCTGCC	GGTCAC	TTAACA	TTATGC	TGAGTG	ATATCC
421	GCGAGC	TCGGTA	CCCGGG	CGAATT	CCAAGC	TT	AGATCT	GGATCC	CCTCTA	GAGTCG
	CGCTCG	AGCCAT	GGGCCC	GCTTAA	GGTTCG	AA	TCTAGA	CCTAGG	GGAGAT	CTCAGC
481	AGGCAT	GCAAGC	TTGGCG	TAATCA	TGGTCA	TAGCTG	TTTCCCT	GTGTGA	AATTGT	TATCCG
	TCCGTA	CGTTTC	AAACCC	ATTAGT	ATCGAC	AAAGGA	CACACT	TTAACA	ATAGCC	
541	CTCACA	ATTCCA	CACAACT	ATACGA	GCCGGA	AGCATA	AAGTGT	AAAGCC	TGGGGT	GCCTAA
	GAGTGT	TAAAGT	GTGTTG	TATGCT	CGGCCCT	TGCTAT	TTTCAG	ATACCA	AGCCCA	CGGATT
601	TGAGTG	AGCTAA	CTCACA	TTAATT	GC GTTG	CGCTCA	CTGCC	GCTTTC	CAGTCG	GGAAAC
	ACTCAC	TCGATT	GAGTGT	AATTAA	CGCAAC	GCGAGT	GACGGG	CGAAAG	GTCAGC	CCTTTG

661	CTGTTC	TGCCAG	CTGCAT	TAATGA	ATCGGC	CAACGC	GCGGGG	AGAGGC	GGTTTG	CGTATT
	GACAGC	ACGGTC	GACGTA	ATTACT	TAGCCG	GTCGCG	GCCTCC	TCTCCG	CCAAAC	GCATAA
721	GGGGCG	TCTTCC	GCTTCC	TCCGTC	ACTGAC	TGCGTG	CGCTCG	GTCGTT	CGGCTG	CGGCGA
	CCCCCG	AGAAGG	CGAAGG	AGCGAG	TGACTG	AGCGAC	GCGAGC	CAGCAA	GCCGAC	GCCGCT
781	GCGGTA	TCAGCT	CACTCA	AAGGCG	GTAATA	CGGTTA	TCCACA	GAATCA	GGGGAT	AACGCA
	CGCCAT	AGTCGA	GTGAGT	TTCCCG	CATTAT	GCCAA	AGGTGT	CTTAGT	CCCCTA	TTGCGT
841	GGAAAG	AACATG	TGAGCA	AAAGGC	CAGCAA	AAGGCC	AGGAAC	CGTAAA	AAGGCC	CGGTTG
	CCTTTC	TTGTAC	ACTCGT	TTTCCG	GTCGTT	TTCCCG	TCCTTG	GCATTT	TGCCGG	GCACAAC
901	CTGGCG	TTTTTC	CATAGG	CTCCGC	CCCCCT	GACGAG	CATCAC	AAAAAT	CGACGC	TCAAGT
	GACCGC	AAAAAG	GTATCC	GAGGCG	GGGGGA	CTGCTC	GTAGTG	TTTTTA	GCTGCG	AGTTCA
961	CAGAGG	TGGCGA	AACCCG	ACAGGA	CTATAA	AGATAA	CAGGCG	TTTCCC	CCTGGA	AGCTCC
	GTCTCC	ACCGTC	TTGGGC	TGTCCT	GATATT	TCTATG	GTCCCG	AAAGGG	GCATTT	GTCCAG
1021	CTCGTG	CGCTCT	CCTGTT	CCGACC	CTGCCG	CTTACC	GGATAC	CTGTCC	GCCTTT	CTCCCT
	GAGCAC	GCGAGA	GGACAA	GGCTGG	GACGGC	GAATGG	CCTATG	GACAGG	CGGAAA	GAGGGA
1081	TCCGGG	AGCGTG	GCGCTT	TCTCAT	AGTCFA	CGCTGT	AGGTAT	CTCAGT	TCGGTG	TAGGTC
	AGCCCT	TCCGAC	CGCGAA	AGAGTA	TCGAGT	GCGACA	TCCATA	GAGTCA	AGCCAC	ATCCAG
1141	GTTTCG	TCCAAG	CTGGGC	TGTGTG	CACGAA	CCCCCC	GTTTCAG	CCCGAC	CGCTGC	CGCTTA
	CAAGCG	AGGTTT	GACCCG	ACACAC	GTGCTT	GGGGGG	CAAGTC	GGGCTG	GCGACG	CGGAAT
1201	TCCGGT	AACTAT	CGTCTT	GAGTCC	AACCCG	GTAAGA	CACGAC	TTATCG	CCACTG	GCAGCA
	AGGCCA	TTGATA	GCAGAA	CTCAGG	TTGGGG	CATTCT	GTGCTG	AAATAGC	GGTGAC	CGTCGT
1261	GCCACT	GGTAAC	AGGATT	AGCAGA	GCGGAG	TATGTA	GGCGGT	GCTACA	GAGTTC	TTGAAG
	CGGTCA	CCATTG	TCCATA	TCTGCT	CGCTCC	ATACAT	CCGCCA	CGATGT	GCACAAG	AACCTC
1321	TGGTGG	CCTAAC	TACGGC	TACACT	AGAAGG	ACAGTA	TTTGGT	ATCTGC	GCTCTG	CTGAAG
	ACCACC	GGATTG	ATGCCG	ATGTGA	TCTTCC	TGTCAT	AAACCA	TAGACG	CGAGAC	GACTTC

1381 CCAGTT ACCTTC GGAAAA AGAGTT GGTAGC TCTTGA TCCGGC AAACAA ACCACC GCTGGT
GGTCAA TGGAAG CCTTTT TCTCAA CCATCG AGAACT AGGCCG TTTGTT TGGTGG CGACCA

1441 AGCGGT GGTFTT TTTGTT FGCAAG CAGCAG ATTACG CGCAGA AAAAAA GGATCT CAAGAA
TCGCCA CCAAAA AAACAA ACCTTC GTCGTC TAAATG GCGTCT TTTTTT CCTAGA GTTCTT

1501 GATCCT TTGATC TTTTCT ACGGGG TCTGAC GCTCAG TGGAAC GAAAAC TCACGT TAAGGG
CTAGGA AACTAG AAAAGA TGCCCC AGACTG CGAGTC ACCTTG CTTTTG AGTGCA ATTCCC

1561 ATTTTG GTCATG AGATTA TCAAAA AGGATC TTCACC TAGATC CTTTTA AATTAA AAATGA
TAAABC CAGTAC TCTAAT AGTTTT TCCTAG AAGTGG ATCTAG GAAAAA TTAAT TTTACT

1621 AGTTTT AAATCA ATCTAA AGTATA TATGAG TAAACT TGGTCT GACAGT TACCAA TGCTTA
TCAAAA TTTAGT TAGATT TCATAT ATACTC ATTTGA ACCAGA CTGTCA ATGGTT ACGAAT

1681 ATCAGT GAGGCA CCTATC TCAGCG ATCTGT CTATTT CGTTCA TCCATA GTTGCC TGAUTC
TAGTCA CTCCGT GGATAG AGTCCG TAGACA GATAAA GCAAGT AGGTAT CAACGG ACTGAG

1741 CCCGTC GTGTAG ATAAC TACGATA CGGGAG GGCTTA CCATCT GGCCCC AGTGCT GCAATG
GGGCAG CACATC TATTGA TGCTAT GCCCTC CCGAAT GGTAGA CCGGGG TCACCA CGTTAC

1801 ATACCG CGAGAC CCACGC TCACCG GCTCCA GATTTA TCAGCA ATAAAC CAGCCA GCCGGA
TATGCC GCTCTG AGTGCC AGTGGT CTAAAT AGTCTG TATTTG GTCGGT GCGCCT

1861 AGGGCC GAGCGC AGAAGT GGTCTT GCAACT TTATCC GCCTCC ATCCAG TCTATT AATGTG
TCCCGG CTCCGC TCTTCA CCAGGA CGTTGA AATAGG CCGGAG TAGGTC AGATAA TTAACA

1921 TGCCCG GAAGCT AGAGTA AGTAGT TCGCCA GTTAAT AGTTTG CGCAAC GTTGTG GCCATT
ACGGCC CTTTCA TCTCAT TCATCA AGCGGT TCAAC CAACAA CGGTAA

1981 GCTACA GGCATC GTGGTG TCACGC TCGTCG TTTGGT ATGGCT TCATTC AGCTCC GGTGCC
CGATGT CCGTAG CACCAC AGTGCC AGCAGC AAACCA TACCGA AGTAAG TCGAGG CCAAGG

2041 CAACGA TCAAGG CGAGTT ACATGA TCCCCC ATGTTG TGCAAA AAAGCG GTTAGC TCCTTC
GTTGCT AGTTCC GCTCAA TGTACT AGGGGG TACAAC ACGTTT TTTCCG CAATCG AGGAAG

14

2101 GGTCTT CCGATC GTTGTG AGAAGT AAGTTG GCCGCA GTGTTA TCACTC ATGGTT ATGGCA
CCAGGA GGCTAG CAACAG TCTTCA TTCAAC CGGCCG CACAAT AGTGAG TACCAA TACCGT

2161 GCAC TG CATAAT TCTCTT ACTGTC ATGCCA TCCGTA AGATGC TTTTCT GTGACT GGTGAG
CGTGAC GTATTA AGAGAA TGACAG TACGGT AGGCAT TCTACG AAAAGA CACTGA CCACTC

2221 TACTCA ACCAAG TCATTC TGAGAA TAGTGT ATGCGG CGACCG AGTTGC TCTTGC CCGGCG
ATGAGT TGGTTC AGTAAG ACTCTT ATCACA TACGCC GCTGGC TCAACG AGAACG GGCCGC

2281 TCAATA CGGGAT AATACC GCGCCA CATAGC AGAACT TAAAAA GTGCTC ATCAT TGGAAA
AGTTAT GCCCTA TTATGG CGCGGT GTATCG TCTTGA AATTTT CACGAG TAGTAA CCTTTT

2341 CGTTCT TCGGGG CGAAAA CTCTCA AGGATC TTACCG CTGTTG AGATCC AGTTCC ATGTAA
GCAAGA AGCCCC GCTTTT GAGAGT TCCTAG AATGGC GACAAC TCTAG TCAAGC TACATT

2401 CCCACT CGTGCA CCCAAC TGATCT TCAGCA TCTTTT ACTTTC ACCAGC GTTTCT GGGTGA
GGGTGA GCACGT GGGTTG ACTAGA AGTCGT AGAAAA TGAAAG TGGTCG CAAAGA CCCACT

2461 GCAAAA ACAGGA AGGCAA AATGCC GCAAAA AAGGGA ATAAGG GCGACA CGGAAA TGTTGA
CGTTTT TGTCCT TCCGTT TTACGG CGTTTT TFCCCT TATTCC CGCTGT GCCTTT ACAACT

2521 ATACTC ATACTC TTCCTT TTTCAA TATTAT TGAAGC ATTTAT CAGGGT TATTGT CTCATG
TATGAG TATGAG AAGGAA AAAGTT ATAATA ACTTCC TAAATA ATAACA GAGTAC

2581 AGCGGA TACATA TTTGAA TGTATT TAGAAA AATAAA CAAATA GGGGT CCGCGC ACATTT
TCGCCT ATGTAT AAACCT ACATAA ATCTTT TTATTT GTTTAT CCCCAG GCGCGC TGTAAA

2641 CCCCAG AAAGTG CCACCT GACGTC TAAGAA ACCATT ATTATC ATGACA TTAACC TATAAA
GGGGCT TTTTCA AAGTGA CTGCGC ATTCCT TTTTAT TTTTAT CCCCAG GCGCGC TGTAAA

2701 AATAGG CGTATC ACGAGG CCCTTT CGTC
TTATCC GCATAG TGCTCC GGGAAA GCAG

15


Protocol for Ligation Using RBC T&A Cloning Vector

1. Centrifuge RBC T&A Cloning Vector and/or PCR DNA tubes to collect contents at the bottom of the tubes.
2. Vortex the ligation buffer vigorously before use.
3. Set up the following items as described below:

Component	User Sample	Positive Control
10X Ligation Buffer A	1 μ l	1 μ l
10X Ligation Buffer B	1 μ l	1 μ l
T&A Cloning Vector	2 μ l	2 μ l
PCR Product	X μ l	***
Control Insert DNA	****	3 μ l
T4 DNA Ligase	1 μ l	1 μ l
Deionized Water	to 10 μ l	to 10 μ l

16

4. Mix the reactions by pipetting.
5. Incubate the reactions for 5 to 15 minutes at 22°C. And, if the maximum of transformants is required, incubate the reactions overnight at 4°C.



17

Notes

1. For blue-white selection, plating should be done onto plates containing 50 µg/ml ampicillin plus 100 µl 100 mM IPTG and 20 µl 50 mg/ml X-gal. Spread the mixture evenly onto a plate at least 2 hours prior to plating transformed cells. After colony growth, the plates can be kept at 4°C and the blue color will continue to develop.
2. Using Colony PCR, clones can be precisely and easily screened.
3. Multiple freeze-thawing DOES NOT affect the quality of RBC T&A Cloning Vector, but exposure to frequent temperature changes will degrade ATP activity in ligase buffer. Always making single-use aliquots of ligase buffers.
4. Pfu DNA polymerase possesses proofreading activity, it does not have the terminal transferase-like activity demonstrated by Taq DNA Polymerase. Ligation reactions using non-tailed amplified DNA will result in no positive colonies. When using non-tailed amplified DNA, refer to methods indicated in note 5 for increasing the ligation efficiency

18

5. Methods for increasing the ligation efficiency:

A. A-tailing	Purified PCR product	X µl
	10X PCR buffer	10 µl
	10 mM dATP	2 µl
	Taq	1 µl

i. Add deionized water to a final volume of 100 µl.

ii. 72°C for 3 hours.

iii. Purify the A-tailed DNA and use in the ligation reaction.

- B. If the maximum of transformants is required, incubate the reactions overnight at 4°C.
- C. The optimized efficiency is using a 1:3 molar ratio of vector DNA to the insert DNA.
- D. Use higher efficiency competent cells e.g. HIT (>10⁸ cfu/µg DNA) series.

19

Troubleshooting

Problem	No Colonies	Less Than 100 Colonies	High Colony Count But High % of Blue with Insert DNA.	No Blue Colonies Are Present on the Plate
Comments	1. Problems in transformation. 2. Low efficiency competent cells.	1. Improper molar ratio of vector DNA to the insert DNA. 2. Bad A-tailed insert DNA.	1. Improper ligation reaction. 2. DNA is inserted, but it's not disrupting the LacZ gene.	1. Ampicillin is inactive. 2. IPTG/X-Gal is insufficient.
Suggestions	Background vector only should yield < 100 colonies, and with A-tailed insert DNA should yield > 500 colonies. If there are no colonies observed on the plate, use high efficiency competent cells to overcome the problem*. *RBC HIT Non-Heat Shock Competent Cells: 1. HIT Competent Cells™-DH5a. 2. HIT Competent Cells™-JM109. 3. HIT Competent Cells™-Blue. 4. HIT Competent Cells™-BL21.	1. Check the quality of A-tailed insert DNA. 2. Use a 1:3 to 1:6 molar ratio of vector DNA to the insert DNA.	1. PCR products should be gel-purified. 2. Avoid multiple freeze-thaw and exposure to frequent temperature changes by making single-use aliquots of Ligase Buffer. 3. If the maximum of transformants is required, incubate the ligation reactions overnight at 4°C. 4. Multiple PCR product, target DNA should be gel-purified before ligation. 5. Use colony PCR to screen the clones.	1. Make sure ampicillin plates are made properly and used within 1 month. 2. Make sure IPTG/X-Gal are fresh. 3. Make sure strain is with blue-screen genotype.

20

21