

RealPfu™ DNA Polymerase

Most Accurate Amplification



Cat. No. RT004

500 units, with dNTP

RealPfu™ DNA Polymerase (5 U/μl): 100 μl

10X Reaction Buffer (with 20 mM Mg²⁺): 2 ml

10 mM dNTPs Mix: 200 μl

Cat. No. RT044

500 units, without dNTP

RealPfu™ DNA Polymerase (5 U/μl): 100 μl

10X Reaction Buffer (with 20 mM Mg²⁺): 2 ml

Recombinant	✓
5' to 3' Exonuclease	✓
3' to 5' Exonuclease	✓
Endonuclease Free	✓

Description

RealPfu™ DNA Polymerase is purified from *Pyrococcus furiosus*. The thermostable enzyme catalyzes the incorporation of nucleotides into duplex DNA in the 5'–3' direction in the presence of Mg²⁺ at 70–80°C. RealPfu™ DNA Polymerase exhibits both 3'–5' exonuclease (proofreading) activity and 5'–3' exonuclease activity, which leads to the lowest error rate and the most reliable high-fidelity PCR.

Unit Definition

One unit of enzyme catalyzes the incorporation of 10 nmol of dNTP into acid-insoluble form in 30 minutes at 74°C.

Error Rate

The error rate of RealPfu™ DNA Polymerase is 1x10⁻⁶ errors per nucleotide per cycle.

Storage Buffer

50 mM Tris-HCl (pH 9.0), 100 mM NaCl, 0.1 mM EDTA, 1% Triton X-100, 5 mM DTT, 50% Glycerol, Stabilizers

10X Reaction Buffer

100 mM KCl, 20 mM MgSO₄ · 7H₂O, 200 mM Tris-HCl (pH 8.8), 1% Triton X-100, 100 mM (NH₄)₂SO₄, 1 mg/ml BSA.

Quality Control

Nuclease activity was not detected after incubation of 1 μg lambda/HindIII DNA with 5 units of RealPfu™ DNA Polymerase in 50 μl reaction volume reaction buffer for 18 hours at 37°C.

General Reaction Conditions

The optimal conditions for the concentration of RealPfu™ DNA Polymerase, primers and template DNA will depend on the system being utilized. It may be necessary to determine the optimal conditions for each individual component.

1. Add the following components to a sterile microtube on ice:

Components	Volume	Final Concentration
10X Reaction Buffer	5 μl	1X
10 mM dNTP Mix	0.5 μl	0.1 μm
Primer Mix (10 μM each)	1 μl	0.2 μm
Template DNA	0.5–10 μl	n/a
RealPfu™ DNA Polymerase (5U/μl)	0.5 μl	2.5 units
ddH ₂ O	add to 50 μl	n/a

2. Suggested Reaction Parameters for RealPfu™ DNA Polymerase:

Segment	Number of Cycles	Temperature	Duration
1	1	94°C	5 minutes
2	25~35	94°C (Denature)	30 seconds
		Primer T _m -5°C (Anneal)	30 seconds
		68°C (Extend)	2 minutes/Kb
3	1	68°C	7 minutes
		4°C	∞

3. Analyze the amplification products by agarose gel electrophoresis and visualize by ethidium bromide staining.

Note: For research use only. Not for use in diagnostic or therapeutic procedures.

