



Real Biotech Corporation

13F-2, No.33, Sec. 1, Minsheng Rd., Banqiao City, Taipei County 220, Taiwan, R. O. C.
Tel: +886 2 2950 9000 Fax: +886 2 2950 0505

PairFast™ Real-Time PCR Mastermix (For Probe System w/ ROX)

Description

PairFast™ Real-Time PCR Mastermix is especially designed for amplifying DNA template below 300 bp in a fast PCR mode. It's supplied as 2-fold concentrated, ready-to-use mixture which is highly sensitive and optimized for use with any real-time PCR cycler and any sequence-specific probes (ex: Taqman, Molecular Beacon...etc). It contains all the factors needed to perform real-time PCR, including RBC SuperiorTaq® HotStart DNA Polymerase, RBC Taqman real-time PCR buffer, ROX passive reference dye, dNTPs and 5mM MgCl₂. The only step to perform real-time PCR is to add the primers, probe and template. With unique composition of RBC SuperiorTaq® HotStart DNA Polymerase and RBC Taqman real-time PCR buffer, PairFast™ Real-Time PCR Mastermix provides fast, highly sensitive detection and accurate quantification. It makes real-time PCR quick, simple and easy.

Specifications

Cat. No.	Product Name	Specification
RT801	PairFast™ Real-Time PCR Mastermix (For Probe System w/ ROX), 100 reactions	2X PairFast™ Real-Time PCR Mastermix: 1ml
RT802	PairFast™ Real-Time PCR Mastermix (For Probe System w/ ROX), 400 reactions	2X PairFast™ Real-Time PCR Mastermix: 4ml
RT803	PairFast™ Real-Time PCR Mastermix (For Probe System w/ ROX), 800 reactions	2X PairFast™ Real-Time PCR Mastermix: 8ml

Features

- High PCR specificity with unique composition of RBC SuperiorTaq® HotStart DNA Polymerase
- Fast, highly sensitive detection and accurate quantification for DNA template below 300 bp
- Optimized, ready-to-use mixture format makes real-time PCR quick, simple and easy

Content

- RBC SuperiorTaq® HotStart DNA Polymerase
- RBC Taqman real-time PCR buffer
- ROX passive reference dye
- dNTP mix including dATP、dCTP、dGTP、dTTP
- 5mM MgCl₂



Real Biotech Corporation

13F-2, No.33, Sec. 1, Minsheng Rd., Banqiao City, Taipei County 220, Taiwan, R. O. C.
Tel: +886 2 2950 9000 Fax: +886 2 2950 0505

Quality Control

Tenfold serial dilution (10^9 ~ 10^{10}) of plasmid DNA were amplified using primers specific to the NNV gene. Triplicate reactions at each concentration were amplified along with no-template controls. Standard curve is $r=0.999$, efficiency=92.4% and standard deviation of $Ct < 1.0$.

Applications

PairFast™ Real-Time PCR Mastermix is compatible with all available real-time cyclers, including instruments from Applied Biosystems, Bio-Rad, Roche...etc. It is optimized for quantitative real-time PCR and two-step RT-PCR using probe detection format.

Shipping and Storage Conditions

PairFast™ Real-Time PCR Mastermix is shipped on dry ice and should be stored immediately upon receipt at -20°C in a constant temperature freezer and protected from light. Avoid repeated freeze–thaw cycles. With proper storage and handled correctly, PairFast™ Real-Time PCR Mastermix can be stored for up to 12 months without showing any deduction in performance and quality.

Protocol

General Reaction Conditions

Our protocol is for a reaction size of 20 μ l. This protocol serves only as a guideline for real-time PCR amplification. Optional reaction conditions may vary and must be individual determined.

Notes: Use disposable tips containing hydrophobic filters to minimize cross-contamination.

1. Prepare the reaction mixture on ice.
2. Thaw the reagents completely, vortex well and then add following components to a sterile microtube on ice:

Component	Volume/ Reaction	Final Concentration
2X PairFast™ Real-Time PCR Mastermix	10 μ l	1X
Forward Primer (10 μ M)	0.6~1.2 μ l	0.3~0.6 μ M
Reverse Primer (10 μ M)	0.6~1.2 μ l	0.3~0.6 μ M
Probe (10 μ M)	0.4~0.8 μ l	0.2~0.4 μ M
RNase-Free Water	Add to 18.0 μ l	

3. Mix above components thoroughly by pipetting up and down and dispense the 18 μ l of mixture into PCR tubes or plates.
4. Add 2 μ l of the DNA or cDNA and mix carefully by pipetting up and down.
5. Suggested Real-Time Cycler Conditions.

Segment	Number of Cycles	Temperature	Time
1	1	95 $^{\circ}$ C	20 seconds
2	40~45	95 $^{\circ}$ C	3 seconds
		58~65 $^{\circ}$ C *	\geq 20 seconds
3	1	4 $^{\circ}$ C	∞

*Optimal annealing temperature is depending on user's primer sequences. Suggested annealing temperature is above the T_m of Primer dimmers, but approximately 3 $^{\circ}$ C below the T_m of the specific PCR product. (T_m dimer < annealing temp. < T_m product).

6. Place the PCR tubes or PCR plates in the thermal cycle and start the cycling program.
7. Perform data analysis of the PCR products.