



Real Biotech Corporation

13F.-2, No.33, Sec. 1, Minsheng Rd., Banqiao City, Taipei County 220, Taiwan, R. O. C.
Tel: +886 2 2950 9000 Fax: +886 2 2950 0505

www.real-biotech.com

HiYield™ Genomic DNA Extraction Kit (Soil)

Cat. No.:	YGSL50	YGSL100
Product Name:	HiYield™ Genomic DNA Extraction Kit (Soil)	
Reactions:	50	100
Sample:	250-500 mg of Soil	
Yield:	Up to 5 µg of pure genomic DNA from 250-500 mg of soil samples	
Elution Volume:	30-100 µl	
Format:	Spin Column	
Operation:	Centrifuge	
Operation Time:	Within 40 Minutes	

Description

HiYield™ Genomic DNA Extraction Kit (Soil) provides a fast and economical method for purification of genomic DNA from microorganisms such as bacteria, archaea, fungi, and algae in soil samples. The soil sample is homogenized using a lysis buffer combined with ceramic beads. Insoluble particles, proteins and PCR inhibitors such as humic acid are then precipitated using a unique inhibitor removal buffer. Residual PCR inhibitors remaining in the clear supernatant are further removed by passing through a specialized PCR inhibitor removal column. Genomic DNA in the sample is then bound by the GD column followed by wash and elution. DNA phenol extraction or alcohol precipitation is not required and the entire procedure can be completed within 40 minutes. The purified genomic DNA is ready for use in PCR, restriction enzyme digestion, and sequencing reactions.

Features

- High-quality genomic DNA from soil sample within 40 minutes.
- Consistently high yields (up to 5 µg of pure genomic DNA from 250-500 mg of soil).
- Complete removal of all contaminants for sensitive downstream applications.
- No phenol, chloroform or alcohol.
- Rapid and simple procedure.

Quality Control

The quality of HiYield™ Genomic DNA Extraction Kit (Soil) is tested on a lot-to-lot basis by isolation of genomic DNA from 250 mg of soil samples. The purified DNA is quantified with a spectrophotometer and the yield of genomic DNA is more than 2 µg with A260/A280 ratio 1.7-2.0. The purified DNA is checked by electrophoresis.



Real Biotech Corporation

13F.-2, No.33, Sec. 1, Minsheng Rd., Banqiao City, Taipei County 220, Taiwan, R. O. C.
Tel: +886 2 2950 9000 Fax: +886 2 2950 0505

www.real-biotech.com

Applications

Purified DNA is ready for direct use in PCR, Southern Blotting, Real-Time PCR, AFLP, RFLP, PADP.

Contents

Item	YGSL50	YGSL100
GSL1 Buffer*	50 ml	85 ml
GSL2 Buffer	15 ml	30 ml
GSL3 Buffer	90 ml	160 ml
Wash Buffer (concentrated)** (Add ethanol)	25 ml (100 ml)	25 ml (100 ml)
Elution Buffer	6 ml	30 ml
Inhibitor Removal Column	50 pcs	100 pcs
GD Column	50 pcs	100 pcs
Bead Tube	50 pcs	100 pcs
2 ml Centrifuge Tube	50 pcs	100 pcs
2 ml Collection Tube	50 pcs	100 pcs

* If precipitates have formed in GSL1 Buffer, warm the buffer in a 37°C water bath, followed by gentle shaking to dissolve.

**Add absolute ethanol (see the bottle label for volume) to Wash Buffer then mix by shaking for a few seconds. Check the box on the bottle showing absolute ethanol is added. Be sure and close the bottle tightly after each use to avoid ethanol evaporation.

Safety Notes

During the procedure, always wear a lab coat, disposable gloves, and protective goggles.

Storage

HiYield™ Genomic DNA Extraction Kit (Soil) shall be shipped and stored dry at room temperature (15-25°C).

Protocol

Please read the entire instruction manual prior to starting.

Things to Do before Starting:

1. If precipitates have formed in GSL1 Buffer, warm the buffer in a 37°C water bath, followed by gentle shaking to dissolve.
2. Add absolute ethanol (see the bottle label for volume) to Wash Buffer then mix by shaking for a few seconds. Check the box on the bottle showing absolute ethanol is added. Be sure and close the bottle tightly after each use to avoid ethanol evaporation.

Additional Requirements:

1.5 ml microcentrifuge tubes, standard vortexer, absolute ethanol.

Caution: During operation, always wear a lab coat, disposable gloves and protective goggles.

Step 1 Sample Lysis	<p>Transfer 250-500 mg of soil to a Bead Tube containing ceramic beads. Add 750 µl of GSL1 Buffer then vortex briefly.</p> <p>NOTE:</p> <ol style="list-style-type: none"> 1. Very dry soil samples can soak up large amounts of GSL1 Buffer. In this case, either reduce the soil amount or add additional GSL1 Buffer to the Bead Tube. For wet soil samples, after transferring to a Bead Tube, centrifuge at 8,000 x g for 1 minute. Remove as much liquid as possible with a pipette before adding GSL1 Buffer. Attach the Bead Tubes horizontally to a standard vortex by taping or using an adapter. Vortex at maximum speed for 10 minutes at room temperature. Centrifuge the Bead Tubes at 8,000 x g for 2 minutes at room temperature to eliminate the foam caused by detergents present in GSL1 Buffer. 2. Preheat the required Elution Buffer (100 µl per sample) to 60°C for DNA elution.
Step 2 PCR Inhibitors Removal	<p>Add 150 µl of GSL2 Buffer to the Bead Tube and vortex for 5 seconds. Incubate at 0-4°C for 5 minutes. Centrifuge at 8,000 x g for 1 minute at room temperature to precipitate insoluble particles and PCR inhibitors. Place an Inhibitor Removal Column (purple ring) in a 2 ml Centrifuge Tube. Transfer 500-600 µl of clear supernatant from the Bead Tube to the Inhibitor Removal Column. Centrifuge at 16,000 x g for 1 minute at room temperature then discard the column. Save the flow-through in the 2 ml Centrifuge Tube for DNA Binding.</p> <p>NOTE: If a pellet is in the flow-through, transfer the clear supernatant to a new 1.5 ml microcentrifuge tube (not provided).</p>

<p>Step 3 DNA Binding</p>	<ol style="list-style-type: none"> 1. Add 900 μl of GSL3 Buffer to the flow-through then mix IMMEDIATELY by shaking vigorously for 5 seconds. 2. Place a GD Column (green ring) in a 2 ml Collection Tube. Transfer 750 μl of sample mixture to the GD Column, centrifuge at 16,000 x g for 1 minute at room temperature then discard the flow-through. Place the GD Column back in the 2 ml Collection Tube. 3. Transfer the remaining sample mixture to the GD Column and centrifuge at 16,000 x g for 1 minute at room temperature. Discard the flow-through then place the GD Column back in the 2 ml Collection Tube.
<p>Step 4 Wash</p>	<ol style="list-style-type: none"> 1. Add 400 μl of GSL3 Buffer to the GD Column. Centrifuge at 16,000 x g for 30 seconds at room temperature. Discard the flow-through then place the GD Column back in the 2 ml Collection Tube. 2. Add 600 μl of Wash Buffer (make sure absolute ethanol was added) to the GD Column. Centrifuge at 16,000 x g for 30 seconds at room temperature. Discard the flow-through then place the GD Column back in the 2 ml Collection Tube. 3. Add 600 μl of Wash Buffer (make sure absolute ethanol was added) to the GD Column and centrifuge at 16,000 x g for 30 seconds at room temperature again. Discard the flow-through then place the GD Column back in the 2 ml Collection Tube. 4. Centrifuge at 16,000 x g for 3 minutes at room temperature to dry the column matrix.
<p>Step 5 DNA Elution</p>	<ol style="list-style-type: none"> 1. Transfer the dried GD Column to a new 1.5 ml microcentrifuge tube. 2. Add 30-100 μl of pre-heated Elution Buffer¹, TE Buffer² or water³ into the CENTER of the column matrix. Let stand for at least 3 minutes to allow Elution Buffer, TE Buffer or water to be completely absorbed. 3. Centrifuge at 16,000 x g for 2 minutes at room temperature to elute the purified DNA. <p>NOTE:</p> <p>¹ If a higher DNA concentration is required, use 30 μl of Elution Buffer (10 mM Tris-HCl, pH8.5 at 25°C) then repeat the Elution step by adding the same 30 μl of Elution Buffer (which now contains the eluted DNA) to the center of the column matrix again. If maximum DNA yield is required, use 100 μl of Elution Buffer.</p> <p>² Using TE (10 mM Tris-HCl, 1 mM EDTA, pH8.0) for elution is beneficial as EDTA preserves DNA for long term storage. However, EDTA will affect PCR and other sensitive downstream applications.</p> <p>³ If using water for elution, ensure the water pH is between 7.0 and 8.5. ddH₂O should be fresh as ambient CO₂ can quickly cause acidification. DNA Eluted in water should be stored at -20°C to avoid degradation.</p>

Troubleshooting

Problem	Possible Reasons / Solution
Low Nucleic Acid Yield	<p><u>Too much starting materials</u></p> <p>Too much soil sample was added to the Bead Tube. Too little space in the Bead Tube results in the beads can not efficiently disrupt the sample.</p>
	<p><u>Sample lysis or homogenization was incomplete</u></p> <p>Horizontally vortex the Bead Tube at the maximum speed using a vortexer at room temperature for 10 minutes or using a Disruptor Genie or similar.</p>
	<p><u>Incorrect DNA elution step</u></p> <p>Pre-heat the Elution Buffer to 60 °C prior to DNA elution. Make sure Elution Buffer is added to the center of the GD Column and is absorbed completely.</p> <p>Ensure to use pre-heated Elution Buffer, TE, or water (60~70°C). If using water for elution, ensure the water pH is between 7.0 and 8.5. ddH₂O should be fresh as ambient CO₂ can quickly cause acidification.</p>
Degraded DNA	<p><u>Inappropriate buffer preparation</u></p> <p>Add appropriate volume of absolute ethanol (see the bottle label for volume) to the Wash Buffer prior to use. If precipitates have formed in GSL1 Buffer, warm the buffer in a 37°C water bath, followed by gentle shaking to dissolve.</p>
	<p><u>Mechanical sample disruption is too strong</u></p> <p>Using alternative lysis method for less DNA shearing: After adding GSL1 Buffer, vortex the Bead Tube at maximum speed for 5 seconds then incubate the Bead Tube at 70 °C for 5 minutes. Repeat these steps for 3 times. This lysis method will reduce DNA shearing but may also reduce DNA yield.</p>
Eluted DNA does not perform well in downstream applications	<p><u>Residual ethanol contamination</u></p> <p>Following the wash step, dry the GD Column with additional centrifugation at 16,000 x g for 5 minutes to ensure the GD Column membrane is completely dry.</p>
	<p><u>PCR inhibitors contamination</u></p> <p>Using diluted DNA (1:10) as a template to reduce the concentration of PCR inhibitors for PCR reactions. The DNA can be further purified using our HiYield Plus™ PCR Purification Kit (cat. no. QPP100) to eliminate PCR inhibitors.</p>