

Cat. No.:	YGT50	YGT100	YGT300
Product Name:	HiYield Genomic DNA Mini Kit (Tissue)		
Reactions:	50	100	300
Sample:	Up to 30 mg of Tissue, Up to 25 mg of Paraffin-Embedded		
	Tissue, Buccal S	Swab and Up to 15 ml	of Amniotic Fluid
Yield:	Up to 50 ug		
Format:	Spin Column		
Operation:	Centrifuge		
<b>Operation Time:</b>		Within 60 Minutes	

# HiYield Genomic DNA Mini Kit (Tissue)

### Introduction

HiYield Genomic DNA Mini Kit (Tissue) is designed specifically for purifying total DNA (including genomic, mitochondrial and viral DNA) from a variety of animal tissue, paraffin-embedded tissue, buccal swab and amniotic fluid. Provided micropestles can efficiently homogenize tissue samples to shorten the time in the Lysis Step. The entire procedure can be completed within 1 hour without phenol/chloroform extraction or alcohol precipitation. The expected yield of genomic DNA is up to 50 µg and the purified DNA (with approximately 20-30 Kb) is suitable for direct use in PCR or other enzymatic reactions.

### Applications

Purified DNA is ready for direct use in PCR, Southern Blotting, Real-Time PCR, AFLP, RFLP, PADP.

#### **Components**

ITEM	YGT50	YGT100	YGT300
GT Buffer	30ml	60ml	155ml
GBT Buffer	40ml	60ml	155ml
W1 Buffer	45ml	45ml	130ml
Wash Buffer*	25ml	25ml	50ml
Elution Buffer	30ml	30ml	75ml
Proteinase K**	11mg	22mg	65mg
GT Column	50pcs	100pcs	300pcs
2 ml Collection Tube	100 pcs	200 pcs	600pcs
Micropestle	50pcs	100pcs	300pcs

\* Add absolute ethanol to the Wash Buffer prior to initial use (see the bottle label for volume).
\*\*Add ddH<sub>2</sub>0 to prepare the Proteinase K (vortex to dissolve and spin down) and store at 4°C (see the bottle label for volume).



### Features

Rapid isolation of ready-to-use DNA. No phenol, chloroform or alcohol. Simple procedure.

#### **Quality Control**

The quality of HiYield Genomic DNA Mini Kit (Tissue) is tested on a lot-to-lot basis by isolation of genomic DNA from 20 mg of mouse liver. Purified DNA is quantified with a spectrophotometer and the yield of genomic DNA is more than 10 ug with A260/A280 ratio 1.7 - 1.9. The purified DNA is checked by electrophoresis.

#### Caution

GBT Buffer contains guanidine hydrochloride which is a harmful irritant. During operation, always wear a lab coat, disposable gloves, and protective goggles.

#### References

(1) Vogelstein, B., and Gillespie, D. (1979) Proc. Natl. Acad. Sci. USA 76, 615.



## **Tissue Protocol**

### Things to do before starting

- 1. Add absolute ethanol to the Wash Buffer prior to initial use (see the bottle label for volume).
- 2. Add ddH<sub>2</sub>0 to prepare the Proteinase K (vortex to dissolve and spin down) and store at 4°C (see the bottle label for volume).

Additional requirements: microcentrifuge tubes, absolute ethanol, RNase A (10 mg/ml), ddH<sub>2</sub>0.

Step 1 Tissue Dissociation	<ul> <li>★Cut up to 30 mg of animal tissue (or 0.5 cm of mouse tail) and transfer it to a 1.5 ml microcentrifuge tube. If the tissue has a higher number of cells (eg. spleen or liver), reduce the starting material to 10 mg.</li> <li>★Use the provided Micropestle to grind the tissue to a pulp.</li> <li>★Add 200 µl of GT Buffer to the tube and continue to homogenize the sample tissue by grinding.</li> </ul>
Step 2 Lysis	<ul> <li>★Add 20 µl of Proteinase K to the sample mixture and mix by vortex.</li> <li>★Incubate at 60°C for 30 minutes to lyse the sample. During incubation, invert the tube every 5 minutes.</li> <li>★Add 200 µl of GBT Buffer and mix by vortex for 5 seconds.</li> <li>★Incubate at 70°C for 20 minutes or until the sample lysate is clear. During incubation, invert the tube every 5 minutes. At this time, preheat the required Elution Buffer (200 µl per sample) to 70°C (for Step 5 DNA Elution).</li> <li>★If there is insoluble material present following incubation, centrifuge for 2 minutes at full speed and transfer the supernatant to a new 1.5 ml microcentrifuge tube.</li> </ul>
	<ul> <li>Optional Step: RNA Degradation</li> <li>If RNA free genomic DNA is required, perform this optional step.</li> <li>★Following 70°C incubation, add 4 µl of RNase A (10 mg/ml) to the sample lysate and mix by vortex.</li> <li>★Incubate at room temperature for 5 minutes</li> </ul>



Step 3 DNA Binding	<ul> <li>★Add 200 µl of absolute ethanol to the sample lysate and vortex immediately for 10 seconds. If precipitate appears, break it up by pipetting.</li> <li>★Place a GT Column in a 2 ml Collection Tube.</li> <li>★Transfer all of the mixture (including any precipitate) to the GT Column.</li> <li>★Centrifuge at full speed for 2 minutes.</li> <li>★Discard the 2 ml Collection Tube containing the flow-through and transfer the GT Column to a new 2 ml Collection Tube.</li> </ul>
Step 4 Wash	<ul> <li>★Add 400 µl of W1 Buffer to the GT Column.</li> <li>★Centrifuge at full speed for 30 seconds.</li> <li>★Discard the flow-through and place the GT Column back in the 2 ml Collection Tube.</li> <li>★Add 600 µl of Wash Buffer (ethanol added) to the GT Column.</li> <li>★Centrifuge at full speed for 30 seconds.</li> <li>★Discard the flow-through and place the GT Column back in the 2 ml Collection Tube.</li> <li>★Centrifuge at for 30 seconds.</li> <li>★Discard the flow-through and place the GT Column back in the 2 ml Collection Tube.</li> <li>★Centrifuge again for 3 minutes at full speed to dry the column matrix.</li> </ul>
Step 5 DNA Elution	<ul> <li>Standard elution volume is 100 µl. If less sample is to be used, reduce the elution volume (30-50 µl) to increase DNA concentration.</li> <li>If higher DNA yield is required, repeat the DNA Elution step to increase DNA recovery and the total elution volume to approximately 200 µl.</li> <li>★Transfer the dried GT Column to a clean 1.5 ml microcentrifuge tube.</li> <li>★Add 100 µl of preheated Elution Buffer or TE to the center of the column matrix.</li> <li>★Let stand for 5 minutes or until the Elution Buffer or TE is absorbed by the matrix.</li> <li>★Centrifuge at full speed for 30 seconds to elute the purified DNA.</li> </ul>



# **Paraffin-Embedded Tissue Protocol**

### Things to do before starting

- 1. Add absolute ethanol to the Wash Buffer prior to initial use (see the bottle label for volume).
- 2. Add ddH<sub>2</sub>0 to prepare the Proteinase K (vortex to dissolve and spin down) and store at 4°C (see the bottle label for volume).

Additional requirements: xylene, microcentrifuge tubes, absolute ethanol, RNase A (10 mg/ml),  $ddH_20$ .

Step 1 Tissue Dissociation	<ul> <li>★Slice small sections (up to 25 mg) from blocks of paraffin-embedded tissue and transfer to a 1.5 ml microcentrifuge tube.</li> <li>★Add 1 ml of xylene to the tube. Vortex vigorously and incubate at room temperature for approximately 10 minutes. Vortex occasionally during incubation.</li> <li>★Centrifuge at full speed for 3 minutes. Remove the supernatant.</li> <li>★Add 1 ml of absolute ethanol to wash the sample pellet and mix by inverting.</li> <li>★Centrifuge at full speed for 3 minutes. Remove the supernatant.</li> <li>★Add 1 ml of absolute ethanol to wash the sample pellet again and mix by inverting.</li> <li>★Centrifuge at full speed for 3 minutes. Remove the supernatant.</li> <li>★Add 1 ml of absolute ethanol to wash the sample pellet again and mix by inverting.</li> <li>★Centrifuge at full speed for 3 minutes. Remove the supernatant.</li> <li>★Add 1 ml of absolute ethanol to wash the sample pellet again and mix by inverting.</li> <li>★Centrifuge at full speed for 3 minutes. Remove the supernatant.</li> <li>★Add 1 ml of absolute ethanol to wash the sample pellet again and mix by inverting.</li> <li>★Centrifuge at full speed for 3 minutes. Remove the supernatant.</li> <li>★Open the tube and Incubate at 37°C for 15 minutes to evaporate any ethanol residue.</li> </ul>
Step 2 Lysis	<ul> <li>★Add 20 µl of Proteinase K to the sample mixture and mix by vortex.</li> <li>★Incubate at 60°C for 30 minutes to lyse the sample. During incubation, invert the tube every 5 minutes.</li> <li>★Add 200 µl of GBT Buffer and mix by vortex for 5 seconds.</li> <li>★Incubate at 70°C for 20 minutes or until the sample lysate is clear. During incubation, invert the tube every 5 minutes. At this time, preheat the required Elution Buffer (200 µl per sample) to 70°C (for Step 5 DNA Elution).</li> <li>★If there is insoluble material present following incubation, centrifuge for 2 minutes at full speed and transfer the supernatant to a new 1.5 ml microcentrifuge tube.</li> </ul>



	<ul> <li>Optional Step: RNA Degradation</li> <li>If RNA free genomic DNA is required, perform this optional step.</li> <li>★Following 70°C incubation, add 4 µl of RNase A (10 mg/ml) to the sample lysate and mix by vortex.</li> <li>★Incubate at room temperature for 5 minutes</li> </ul>
Step 3 DNA Binding	<ul> <li>★Add 200 µl of absolute ethanol to the sample lysate and vortex immediately for 10 seconds. If precipitate appears, break it up by pipetting.</li> <li>★Place a GT Column in a 2 ml Collection Tube.</li> <li>★Transfer all of the mixture (including any precipitate) to the GT Column.</li> <li>★Centrifuge at full speed for 2 minutes.</li> <li>★Discard the 2 ml Collection Tube containing the flow-through and transfer the GT Column to a new 2 ml Collection Tube.</li> </ul>
Step 4 Wash	<ul> <li>★Add 400 µl of W1 Buffer to the GT Column.</li> <li>★Centrifuge at full speed for 30 seconds.</li> <li>★Discard the flow-through and place the GT Column back in the 2 ml Collection Tube.</li> <li>★Add 600 µl of Wash Buffer (ethanol added) to the GT Column.</li> <li>★Centrifuge at full speed for 30 seconds.</li> <li>★Discard the flow-through and place the GT Column back in the 2 ml Collection Tube.</li> <li>★Centrifuge at for 30 seconds.</li> <li>★Discard the flow-through and place the GT Column back in the 2 ml Collection Tube.</li> <li>★Centrifuge again for 3 minutes at full speed to dry the column matrix.</li> </ul>
Step 5 DNA Elution	Standard elution volume is 100 µl. If less sample is to be used, reduce the elution volume (30-50 µl) to increase DNA concentration. If higher DNA yield is required, repeat the DNA Elution step to increase DNA recovery and the total elution volume to approximately 200 µl. *Transfer the dried <b>GT Column</b> to a clean 1.5 ml microcentrifuge tube. *Add <b>100 µl of preheated Elution Buffer</b> or TE to the center of the column matrix. *Let stand for 5 minutes or until the <b>Elution Buffer</b> or TE is absorbed by the matrix. *Centrifuge at full speed for 30 seconds to elute the purified DNA.



## **Buccal Swab Protocol**

### Things to do before starting

- 1. Add absolute ethanol to the Wash Buffer prior to initial use (see the bottle label for volume).
- 2. Add ddH<sub>2</sub>0 to prepare the Proteinase K (vortex to dissolve and spin down) and store at 4°C (see the bottle label for volume).

Additional requirements: Swab (cotton, DACRON or C.E.P. swabs), PBS (phosphate-buffered saline), microcentrifuge tubes, absolute ethanol, ddH<sub>2</sub>0.

Step 1 Sample Preparation	★Scrape the swab firmly against the inside of each cheek 6-7 times and air-dry the swab. (The person providing the sample should not eat or drink for at least 30 minutes prior to sample collection to avoid contamination.)
Step 2 Lysis	<ul> <li>★Add 500 µl of GT Buffer and 20 µl of Proteinase K to a 1.5 ml microcentrifuge tube.</li> <li>★Place the buccal swab into the tube and incubate at 60°C for 10 minutes.</li> <li>★Discard the swab and add 500 µl of GBT Buffer to the lysate.</li> <li>★Vortex immediately and incubate at 60°C for 10 minutes.</li> <li>★At this time, preheat the required Elution Buffer (200 µl per sample) in a 60°C water bath (for Step 5 DNA Elution).</li> </ul>
Step 3 DNA Binding	<ul> <li>★Add 500 µl of absolute ethanol to the sample lysate and vortex immediately.</li> <li>★Place a GT Column in a 2 ml Collection Tube.</li> <li>★Transfer 700 µl of the mixture to the GT Column.</li> <li>★Centrifuge at full speed for 1 minute.</li> <li>★Repeat the DNA Binding Step by transferring the remaining mixture to the GT Column.</li> </ul>



	★Add 400 μI of W1 Buffer to the GT Column.
Stop 4	★Centrifuge at full speed for 30 seconds.
	$\bigstar$ Discard the flow-through and place the <b>GT Column</b> back in the
	2 ml Collection Tube.
Wash	★Add 600 μl of Wash Buffer (ethanol added) to the GT Column.
Theorem 1	★Centrifuge at full speed for 30 seconds.
	★ Discard the flow-through and place the <b>GT Column</b> back in the
	2 ml Collection Tube.
	$\star$ Centrifuge again for 3 minutes at full speed to dry the column matrix.
	Standard elution volume is 100 ul. If less sample is to be used, reduce the
	elution volume (30-50 ul) to increase DNA concentration.
	If higher DNA yield is required, repeat the DNA Elution step to increase
	DNA recovery and the total elution volume to approximately 200 µl.
Step 5	
<b>DNA Elution</b>	$\star$ Transfer the dried <b>GT Column</b> to a clean 1.5 ml microcentrifuge tube.
	★Add <b>100 µl of preheated Elution Buffer</b> or TE to the center of the column matrix
	$\star$ Let stand for 5 minutes or until the <b>Elution Buffer</b> or TE is absorbed by
	the matrix.
	$\star$ Centrifuge at full speed for 30 seconds to elute the purified DNA.



# **Amniotic Fluid Protocol**

### Things to do before starting

- 1. Add absolute ethanol to the Wash Buffer prior to initial use (see the bottle label for volume).
- 2. Add ddH<sub>2</sub>0 to prepare the Proteinase K (vortex to dissolve and spin down) and store at 4°C (see the bottle label for volume).

Additional requirements: centrifuge tubes, microcentrifuge tubes, absolute ethanol, RNase A (10 mg/ml),  $ddH_20$ .

Step 1 Cell Harvesting/ Prelysis	<ul> <li>★Transfer approximately 10 ml (up to 15 ml) of amniotic fluid to a 15 ml centrifuge tube.</li> <li>★Centrifuge for 3 minutes at full speed and discard the supernatant.</li> </ul>
Step 2 Lysis	<ul> <li>★Add 200 µl of GT Buffer to resuspend the pellet and transfer the mixture to a 1.5 ml microcentrifuge tube.</li> <li>★Add 10 µl of Proteinase K (10 mg/ml) to the sample mixture and mix by vortex.</li> <li>★Incubate at 60°C for 30 minutes to lyse the sample. During incubation, invert the tube every 5 minutes.</li> <li>★Add 200 µl of GBT Buffer to the 1.5 ml microcentrifuge tube and mix by vortex for 5 seconds.</li> <li>★Incubate at 70°C for 20 minutes or until the sample lysate is clear. During incubation, invert the tube every 5 minutes.</li> <li>★At this time, preheat the required Elution Buffer (200 µl per sample) to 70°C. (for Step 5 DNA Elution)</li> </ul>
	<ul> <li>Optional Step: RNA Degradation</li> <li>If RNA free genomic DNA is required, perform this optional step.</li> <li>★Following 70°C incubation, add 5 µl of RNase A (10 mg/ml) to the sample lysate and mix by vortex.</li> <li>★Incubate at room temperature for 5 minutes</li> </ul>



Step 3 DNA Binding	<ul> <li>★Add 200 µl of absolute ethanol to the sample lysate and vortex immediately for 10 seconds. If precipitate appears, break it up by pipetting.</li> <li>★Place a GT Column in a 2 ml Collection Tube.</li> <li>★Transfer all of the mixture (including any precipitate) to the GT Column.</li> <li>★Centrifuge at full speed for 1 minute.</li> <li>★Discard the 2 ml Collection Tube containing the flow-through and transfer the GT Column to a new 2 ml Collection Tube.</li> </ul>
Step 4 Wash	<ul> <li>★Add 400 µl of W1 Buffer to the GT Column.</li> <li>★Centrifuge at full speed for 30 seconds.</li> <li>★Discard the flow-through and place the GT Column back in the 2 ml Collection Tube.</li> <li>★Add 600 µl of Wash Buffer (ethanol added) to the GT Column.</li> <li>★Centrifuge at full speed for 30 seconds.</li> <li>★Discard the flow-through and place the GT Column back in the 2 ml Collection Tube.</li> <li>★Centrifuge at flow-through and place the GT Column back in the 2 ml Collection Tube.</li> <li>★Centrifuge again for 3 minutes at full speed to dry the column matrix.</li> </ul>
Step 5 DNA Elution	<ul> <li>Standard elution volume is 100 µl. If less sample is to be used, reduce the elution volume (30-50 µl) to increase DNA concentration.</li> <li>If higher DNA yield is required, repeat the DNA Elution step to increase DNA recovery and the total elution volume to approximately 200 µl.</li> <li>★Transfer the dried GT Column to a clean 1.5 ml microcentrifuge tube.</li> <li>★Add 100 µl of preheated Elution Buffer or TE to the center of the column matrix.</li> <li>★Incubate the GT Column at 37°C in an incubator for 10 minutes.</li> <li>★Centrifuge at full speed for 1 minute to elute the purified DNA.</li> </ul>



# Hairs That Include Follical Cells Protocol

### Things to do before starting

- 1. Add absolute ethanol to the Wash Buffer prior to initial use (see the bottle label for volume).
- 2. Add ddH<sub>2</sub>0 to prepare the Proteinase K (vortex to dissolve and spin down) and store at 4°C (see the bottle label for volume).

Additional requirements: microcentrifuge tubes, absolute ethanol, RNase A (10 mg/ml), ddH<sub>2</sub>0.

Step 1 Sample Lysis	<ul> <li>★Add 10 hairs that include follical cells to a 2 ml Collection Tube.</li> <li>★Add 200 µl of GT Buffer and 20 ul of Proteinase K to the 2 ml Collection Tube and mix by shaking. (be sure the buffer completely cover the sample)</li> <li>★Incubate at 60°C for 30 minutes to lyse the sample. During incubation, shaking the 2 ml Collection Tube every 5 minutes.</li> <li>★Add 200 ul of GBT Buffer and mix vigorously.</li> <li>★Incubate at 70°C for 20 minutes and during incubation, shaking the 2 ml Collection Tube every 5 minutes.</li> <li>★At this time, preheat the required Elution Buffer (100 µl per sample) to 70°C (for DNA Elution step).</li> <li>★After incubation, centrifuge for 5 minutes at 14,000 x g and transfer the supernatant to a new 2 ml Collection Tube.</li> </ul>
	<ul> <li>Optional Step: RNA Degradation</li> <li>If RNA free genomic DNA is required, perform this optional step.</li> <li>★Following 70°C incubation, add 4 µl of RNase A (10 mg/ml) to the sample lysate and mix by vortex.</li> <li>★Incubate at room temperature for 5 minutes</li> </ul>
Step 2 DNA Binding	<ul> <li>★Add 200 µl of absolute ethanol to the sample lysate in the 2 ml Collection Tube from the previous step.</li> <li>★Mix immediately by pipetting 5-10 times.</li> <li>★Place a GT Column on a new 2 ml Collection Tube.</li> <li>★Transfer the lysate mixture to the GT Column.</li> <li>★Centrifuge for 5 minutes at 14,000 x g.</li> </ul>



Step 3 Wash	<ul> <li>★Add 400 µl of W1 Buffer to the GT Column.</li> <li>★Centrifuge for 5 minutes at 14,000 x g.</li> <li>★Add 600 µl of Wash Buffer (ethanol added) to the GT Column to wash again.</li> <li>★Centrifuge for 5 minutes at 14,000 x g.</li> <li>★Discard the flow-through and place the GT Column back on the 2 ml Collection Tube.</li> <li>★Centrifuge for 10 minutes at 14,000 x g to remove any ethanol residue</li> </ul>
Step 4 DNA Elution	<ul> <li>Standard elution volume is 100 µl. If less sample is to be used, reduce the elution volume (50-100 µl) to increase DNA concentration.</li> <li>If higher DNA yield is required, repeat the DNA Elution step to increase DNA recovery and the total elution volume to approximately 200 µl.</li> <li>★ Transfer the GT Column to a clean 1.5 ml microcentrifuge tube</li> <li>★ Add 50 µl of preheated Elution Buffer or TE to the center of the GT Column.</li> <li>★ Let stand for 3 minutes or until the Elution Buffer or TE is absorbed by the matrix.</li> <li>★ Centrifuge for 5 minutes at 14,000 x g to elute the purified DNA.</li> <li>★ Add 50 µl of preheated Elution Buffer or TE to the center of the GT Column.</li> <li>★ Let stand for 3 minutes or until the Elution Buffer or TE is absorbed by the matrix.</li> <li>★ Centrifuge for 5 minutes at 14,000 x g to elute the purified DNA.</li> <li>★ Add 50 µl of preheated Elution Buffer or TE is absorbed by the matrix.</li> <li>★ Centrifuge for 5 minutes at 14,000 x g to elute the purified DNA.</li> <li>★ Let stand for 3 minutes or until the Elution Buffer or TE is absorbed by the matrix.</li> <li>★ Centrifuge for 5 minutes at 14,000 x g to elute the purified DNA.</li> <li>★ Let stand for 3 minutes or until the Elution Buffer or TE is absorbed by the matrix.</li> <li>★ Centrifuge for 5 minutes at 14,000 x g to elute the purified DNA.</li> <li>★ Use 5-10 µl for PCR or qPCR (dependent on how many volume of Elution buffer you use)</li> </ul>



### Troubleshooting

Problem	Possible Reason/ Solution
Column clogged	Too much sample was used. ★Reduce sample volume or separate into multiple tubes.
	<ul> <li>Sample tissue was not lysed completely.</li> <li>★Add additional Proteinase K and extend the incubation time in the Lysis Step.</li> <li>★Following the Lysis Step, centrifuge for 2 minutes at full speed to remove sample debris. Transfer the supernatant to a new microcentrifuge tube and proceed with the DNA Binding Step.</li> </ul>
	<ul> <li>Precipitate was formed at DNA Binding Step.</li> <li>★Reduce the sample material.</li> <li>★Prior to loading the column, break up precipitate in ethanol-added lysate.</li> </ul>
Low yield	<ul> <li>Sample tissue was not lysed completely</li> <li>★Add additional Proteinase K and extend the incubation time in the Lysis Step.</li> </ul>
	<ul> <li>Column was clogged at DNA Binding Step</li> <li>★ Following the Lysis Step, remove the insoluble debris by centrifugation.</li> <li>★ Prior to loading the column, break up the precipitate in the ethanol added lysate.</li> </ul>
	Incorrect DNA Elution Step. ★Ensure that Elution Buffer was added and absorbed to the center of GT Column matrix completely.
	Incomplete DNA Elution. ★Elute twice to increase yield.
Eluted DNA does not perform well in downstream applications.	<ul> <li>Residual ethanol contamination.</li> <li>★ Following the wash step, dry GT Column with additional centrifugation at full speed for 5 minutes or incubation at 60°C for 5 minutes.</li> </ul>
	RNA contamination. ★Perform Optional RNA degradation Step.
	Genomic DNA was degraded. ★Use fresh samples or freeze fresh samples in liquid nitrogen immediately and store at -80°C.