



Real Biotech Corporation

13F.-2, No.33, Sec. 1, Minsheng Rd., Banqiao City, Taipei County 220, Taiwan, R. O. C.
Tel: +886 2 2950 9000 Fax: +886 2 2950 0505

www.real-biotech.com

HiYield Genomic DNA Isolation Kit (Yeast)

| Cat. No.: | YGY100 | YGY300 |
|-----------------|---|--------|
| Product Name: | HiYield Genomic DNA Isolation Kit (Yeast) | |
| Reactions: | 100 | 300 |
| Sample: | Up to 2×10^8 of yeast and a variety of fungus species | |
| Yield: | Up to 14.5ug for 2×10^8 of <i>Saccharomyces cerevisiae</i> | |
| Elution Volume: | 50-100 μ l | |
| Format: | Reagent | |
| Operation: | Centrifuge | |
| Operation Time: | Within 40 Minutes | |

Description

HiYield Genomic DNA Isolation Kit (Yeast) is ideal for purification of high-molecular-weight genomic, mitochondrial or viral DNA from *Saccharomyces cerevisiae* and a variety of other yeast and fungus species in a simple and gentle reagent DNA precipitation method. The convenient, scalable purification system removes contaminants and enzyme inhibitors such as proteins and divalent cations. Purified DNA ($A_{260}/A_{280}=1.8-2.0$) is ready for immediate use in sensitive downstream applications or for archiving.

Features

1. Convenient, scalable purification procedure.
2. Reproducible recoveries for constant results.
3. Complete removal of all contaminants for sensitive downstream applications.

Contents

| Item | YGY100 | YGY300 |
|------------------------|--------|--------|
| Cell Suspension Buffer | 90 ml | 225 ml |
| Cell Lysis Buffer | 40 ml | 100 ml |
| Protein Removal Buffer | 15 ml | 40 ml |
| DNA Hydration Buffer* | 50 ml | 50 ml |

* DNA Hydration Buffer: 10 mM Tris-HCl, 1 mM EDTA, pH 8.0.

Storage

HiYield Genomic DNA Isolation Kit (Yeast) shall be shipped and stored dry at room temperature (15-25°C). With proper storage, HiYield Genomic DNA Isolation Kit (Yeast) can be stored for up to 12 months without showing any deduction in performance and quality.



Applications

Purified DNA is highly stable and ready for use in a wide range of applications, such as: PCR, AFLP, RFLP/PADP, Southern Blotting, Real-Time PCR.

Quality Control

The quality of HiYield Genomic DNA Isolation Kit (Yeast) is tested on a lot-to-lot basis by isolation of Genomic DNA from 2×10^8 of *Saccharomyces cerevisiae*. A 15 μ l aliquot of purified genomic DNA from a 100 μ l eluate is analyzed by electrophoresis on a 1% agarose gel. Genomic DNA is quantified with a spectrophotometer and the yield of Genomic DNA is more than 10 μ g.

Important Notes

Caution:

HiYield Genomic DNA Isolation Kit (Yeast) contains irritants. During operation, always wear a lab coat, disposable gloves, protective goggles and (anti-fog) procedure mask.

Reagents to Be Supplied by User:

1.5 ml microcentrifuge tubes, β -mercaptoethanol, zymolase or lyticase, RNase A (50 mg/ml), isopropanol, absolute ethanol for preparing 70% ethanol in ddH₂O.

DNA Hydration Buffer:

Using DNA Hydration Buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0) is beneficial as EDTA preserves DNA for long term storage. However, EDTA will affect PCR and other sensitive downstream applications. If using water instead of DNA Hydration Buffer, ensure the water pH is between 7.0 and 8.5. ddH₂O should be fresh as ambient CO₂ can quickly cause acidification. DNA in water should be stored at -20°C to avoid degradation.

Yeast Double-stranded RNA Killer Virus:

Yeast Double-stranded RNA Killer Virus can be found in the majority of *S. cerevisiae* laboratory strains. The ds RNA Killer Virus produces a ~5000 bp band, sometimes an extremely faint ~2000 bp band, on an agarose gel. These bands do not interfere with most downstream applications, excluding exact quantitation of the sample. The ds RNA Killer Virus can be cured by decreasing the suggested amount of isopropanol addition by $\frac{1}{2}$ in DNA Precipitation step; however, this will also result in reduced genomic DNA yield.

Yield and Quality of Purified DNA:

HiYield Genomic DNA Isolation Kit (Yeast) is designed to purify high yields of high-quality DNA. The actual yield depends on the sample type, genome size of the source, and the quality of the starting material.

Protocol

Please read the entire important notes prior to starting.

| | |
|---------------------------------------|---|
| Step 1 Cell Harvesting | <ol style="list-style-type: none"> 1. Transfer fungus cells (up to 2×10^8) to a 1.5 ml microcentrifuge tube. 2. Harvest fungus cells by centrifugation for 10 minutes at 5,000 x g. 3. Discard the supernatant and resuspend the pellet in 600 μl of Cell Suspension Buffer. 4. Add 2 μl of β-mercaptoethanol and 200 U of lyticase or zymolase. Then incubate at 30°C for 30 minutes. 5. Centrifuge the mixture for 10 minutes at 2,000 x g to harvest the spheroplast and then remove the supernatant. |
| Step 2 Cell Lysis | <ol style="list-style-type: none"> 1. Add 300 μl of Cell Lysis Buffer and resuspend the cell pellet by pipette. 2. Incubate at 60°C for at least 10 minutes to ensure the sample lysate is clear. During incubation, invert the tube every 3 minutes. <p><u>Optional</u> RNA Removal Step: Following 60°C incubation, add 5 μl of RNase A (50 mg/ml) to the clear sample lysate and mix by vortex. Incubate at room temperature for 10 minutes.</p> |
| Step 3 Protein Removal | <ol style="list-style-type: none"> 1. Add 100 μl of Protein Removal Buffer to the sample lysate and vortex IMMEDIATELY for 10 seconds. 2. Centrifuge at 14,000 -16,000 x g for 3 minutes to form a tight and white protein pellet. If the pellet is not tight enough, incubate on ice for 5 minutes followed by centrifugation at 14,000 -16,000 x g for another 3 minutes. |
| Step 4 DNA Precipitation | <ol style="list-style-type: none"> 1. Being careful not to draw any of the protein pellet into the pipette, transfer the supernatant from Step 3 to a new 1.5 ml microcentrifuge tube. 2. Add 300 μl of isopropanol and mix well by gently inverting 20 times. 3. Centrifuge at 14,000 -16,000 x g for 5 minutes. 4. Carefully remove the supernatant and add 300 μl of 70% ethanol to wash the pellet. 5. Centrifuge at 14,000 -16,000 x g for 3 minutes. 6. Discard the supernatant and then air-dry the pellet for 10 minutes. (DO NOT dry the DNA pellet with by vacuum centrifuge and avoid over drying the DNA pellet.) |
| Step 5 DNA Hydration | <ol style="list-style-type: none"> 1. Add 50-100 μl of DNA Hydration Buffer or ddH₂O and then incubate at 60°C for 10 minutes to dissolve the DNA pellet. (Occasionally tapping the bottom of the tube during incubation will promote DNA rehydration.) |