# www.real-biotech.com

Store at ୵୷ୡଂ୲

## Lyophilized Proteinase K

Molecular Biology Grade

Cat. No. YPK11 Proteinase K: 11 mg Cat. No. YPK65 Proteinase K:65 ma Cat. No. YPK100 Proteinase K: 100 ma Cat. No. YPK500 Proteinase K: 500 mg Cat. No. YPK1G Proteinase K:1g

## Format: Lyophilized powder

Source: Tritirachium album limber (Recombinant) Purity: Over 99.5% (Native-PAGE and laser scan assay) Activity: >30 units/mg protein (hemoglobin, pH7.5, 37°C)

## Unit definition

One unit liberates Folin-positive amino acids and peptides corresponding to 1 umol tyrosine in 1 minute at 37°C using denatured hemoalobin as substrate.

## **Storage Conditions**

Lyopholized Proteinase K shall be stored at 4°C in a constant refrigerator and protected from moisture. Proteinase K solutions are stable at 2-8°C for 6 months. Store lyopholized Proteinase K at -20°C for long-term storage.



#### Description

Proteinase K is a stable and highly reactive serine protease. It is active with SDS, urea and EDTA and it is strongly inhibited by PMSF. Proteinase K remains active in detergents and reagents over a broad pH range (4.0-12.5, optimum pH8.0) and is also stable over the temperature range of 25°C to 65°C during use. Proteinase K is useful for protein digestion and contamination removal during nucleic acid purification to prepare RNA and high molecular weight DNA for subsequent reactions

## **Preparation Instructions:**

- 1. Carefully open the Proteinase K vial.
- 2. Add sterilized ddH<sub>2</sub>O as indicated on the vial label.
- 3. Gently mix to dissolve the Proteinase K powder.
- 4. Store the Proteinase K solution at 2-8°C
- 5. Transfer up to 25 mg of animal tissue to a 1.5 ml microcentrifuge tube.
- 6. Add 200 µl of lysis buffer (e.g. 30 mM Tris-Cl, 10 mM EDTA, 1% SDS) to the sample.
- 7. Add 20 μl of Proteinase K solution to the sample then mix by vortex.
- 8. Incubate the sample at 60°C for 2-3 hours or until the sample lysate becomes clear.

#### Note

Maintaining consistent storage conditions is essential as temperature fluctuations can affect product stability. For optimal performance, prepare Proteinase K just prior to use.

## Applications

Proteinase Kis active over the pH range 4~12. It can be used in: (1) Proteolytic inactivation of nucleases during the isolation of DNA/RNA

(2) Inactivation of RNases, DNases and enzymes in reactions, (3) Improving cloning efficiency of PCR products, (4) Preparation of tissue sections for in situ hybridization...etc.

#### **Quality Control**

DNases / RNases / Exonucleases: (not detected).

#### Related Products

YRI001 RNAarmor<sup>™</sup> RNase Inhibitor (2,000 U) YRD100 RNAarmor™RNase/DNaseRemovalReagent(100ml) YRD500 RNAarmor<sup>™</sup> RNase/DNase Removal Reagent (500 ml) YRE500 RNAarmor™ RNase Removal Reagent (500 ml) YRS100 RNAstill<sup>™</sup> RNA Stabilization Reagent (100 ml) YRS250 RNAstilI<sup>™</sup> RNA Stabilization Reagent (250 ml) YRS500 RNAstill™RNA Stabilization Reagent (500 ml) YRR001 RNAstill™ RNA Storage Solution (1 ml) YRR010 RNAstill™ RNA Storage Solution (10ml)

#### For research use only

Not intended for any animal or human therapeutic or diagnostic use. ©2018 Real Biotech Corporation. All rights reserved.

#### www.real-biotech.com

# Lyophilized Proteinase K Molecular Biology Grade

Cat. No. YPK11 Proteinase K:11 mg Cat. No. YPK65 Proteinase K:65 mg Cat. No. YPK100 Proteinase K: 100 ma

Cat. No. YPK500 Proteinase K: 500 mg Cat. No. YPK1G Proteinase K:1 a

Format: Lyophilized powder Source: Tritirachium album limber (Recombinant) Purity: Over 99.5% (Native-PAGE and laser scan assay)

Activity: >30 units/mg protein (hemoglobin, pH7.5, 37°C)

#### Unit definition

One unit liberates Folin-positive amino acids and peptides corresponding to 1 µmol tyrosine in 1 minute at 37°C using denatured hemoglobin as substrate.

## Storage Conditions

Lyopholized Proteinase K shall be stored at 4°C in a constant refrigerator and protected from moisture. Proteinase K solutions are stable at 2-8°C for 6 months. Store lyopholized Proteinase K at -20°C for long-term storage.

### Description

Proteinase K is a stable and highly reactive serine protease. It is active with SDS, urea and EDTA and it is strongly inhibited by PMSF. Proteinase K remains active in detergents and reagents over a broad pH range (4.0-12.5, optimum pH8.0) and is also stable over the temperature range of 25°C to 65°C during use. Proteinase K is useful for protein digestion and contamination removal during nucleic acid purification to prepare RNA and high molecular weight DNA for subseauent reactions.

#### Preparation Instructions:

- 1. Carefully open the Proteinase K vial.
- 2. Add sterilized ddH2O as indicated on the vial label.
- 3. Gently mix to dissolve the Proteinase K powder.
- 4. Store the Proteinase K solution at 2-8°C.
- 5. Transfer up to 25 mg of animal tissue to a 1.5 ml microcentrifuge tube.
- 6. Add 200 µl of lysis buffer (e.g. 30 mM Tris-Cl, 10 mM EDTA, 1% SDS) to the sample. 7. Add 20 µl of Proteinase K solution to the sample
- then mix by vortex.
- 8. Incubate the sample at 60°C for 2-3 hours or until the sample lysate becomes clear.

#### Note

Maintaining consistent storage conditions is essential as temperature fluctuations can affect product stability. For optimal performance, prepare Proteinase Kjust prior to use.

## **Applications**

Proteinase Kis active over the pH range 4~12. It can be used in: (1) Proteolytic inactivation of nucleases during the isolation of DNA/RNA,

(2) Inactivation of RNases, DNases and enzymes in reactions, (3) Improving cloning efficiency of PCR products, (4) Preparation of tissue sections for in situ hybridization...etc.

#### **Ouality Control**

DNases / RNases / Exonucleases: (not detected).

#### **Related Products**

YRI001 RNAarmor<sup>™</sup> RNase Inhibitor (2,000 U) YRD100 RNAarmor™RNase/DNaseRemovalReagent(100ml) YRD500 RNAarmor™ RNase/DNase Removal Reagent (500ml) YRE500 RNAarmor™ RNase Removal Reagent (500 ml) YRS100 RNAstill™ RNA Stabilization Reagent (100ml) YRS250 RNAstill™ RNA Stabilization Reagent (250ml) YRS500 RNAstill<sup>™</sup>RNA Stabilization Reagent (500 ml) YRR001 RNAstill™ RNA Storage Solution (1 ml) YRR010 RNAstill™ RNA Storage Solution (10ml)



Not intended for any animal or human therapeutic or diagnostic use.

For research use only.

©2018 Real Biotech Corporation. All rights reserved.

